plasma care®

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Dejonckheere, C. S. et al. Non-Invasive Physical Plasma for Preventing Radiation Dermatitis in Breast Cancer: A First-In-Human Feasibility Study. *Pharmaceutics* 14, 1767 (2022).

Dejonckheere, C. S. et al. Non-invasive physical plasma for preventing radiation dermatitis in breast cancer: Results from an inpatient-randomised double-blind placebo-controlled trial. *Clinical and Translational Radiation Oncology* 44, 100699 (2024).

Nakayama, B., Garcia, L. T. & Serena, T. Efficacy of Cold Atmospheric Plasma in Chronic Diabetic Foot Ulcer Management: A Case Report. *Am J Case Rep* 25, (2024).

Hämmerle, G., Ascher, S. & Gebhardt, L. Positive effects of cold atmospheric plasma on pH in wounds: a pilot study. *J Wound Care* 32, 530–536 (2023).

Ligresti, C. Use of Cold Plasma in the Treatment of Infected Wounds. *JSRP* (2024)

Ricci, E., Pittarello, M., Ricci, A., Ascher, S. & Gebhardt, L. Treatment of infected wounds with cold atmospheric plasma: a case series. *J Wound Care* 34, 506–512 (2025).

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Terabe, Y., Kaneko, N. & Ando, H. Treating hard-to-heal skin and nail onychomycosis of diabetic foot with plasma therapy. *Dermatol Ther* 34, e15127 (2021).

Terabe, Y., Kaneko, N., Nakabayashi, K. & Ando, H. Using Plasma to Treat Chronic Foot Ulcer Infections. *International Journal of Surgical Wound Care* 3, 33–36 (2022).

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Neonatal and pediatric wound care A contemporary perspective on innovations and best practices Guido CIPRANDI et. al. - Minerva Medica - Books - Medicine - Pediatrics.

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Eichler, J., Rulik, B., Abazid, A. & Stope, M. B. Non-invasive physical plasma improves conventional wound management of cut and bite wounds in wild European hedgehogs. *Sci Rep* 15, 2744 (2025).

Reichold, L. Z. et al. Cellular Response of Immune Cells in the Upper Respiratory Tract After Treatment with Cold Atmospheric Plasma In Vitro. *Int J Mol Sci* 26, 255 (2024).

Karrer, S. et al. In Vitro Safety Study on the Use of Cold Atmospheric Plasma in the Upper Respiratory Tract. *Cells* 13, 1411 (2024).

Kupke, L. S. et al. Cold Atmospheric Plasma Promotes the Immunoreactivity of Granulocytes In Vitro. *Biomolecules* 11, 902 (2021).

Literatur

- ¹Zimmermann, J. L. et al. Test for bacterial resistance build-up against plasma treatment. *New Journal of Physics* 14, 073037 (2012).
- ²Isbary, G. et al. A first prospective randomized controlled trial to decrease bacterial load using cold atmospheric argon plasma on chronic wounds in patients. *Br. J. Dermatol.* 163, 78–82 (2010).
- ³Daeschlein, G. et al. Skin and wound decontamination of multidrug-resistant bacteria by cold atmospheric plasma coagulation. *J Dtsch Dermatol Ges* 13, 143–150 (2015).
- ⁴Arndt, S. et al. Effects of cold atmospheric plasma (CAP) on β -defensins, inflammatory cytokines, and apoptosis-related molecules in keratinocytes in vitro and in vivo. *PLoS ONE* 10, e0120041 (2015).
- ⁵Arndt, S. et al. Cold atmospheric plasma (CAP) changes gene expression of key molecules of the wound healing machinery and improves wound healing in vitro and in vivo. *PLoS ONE* 8, e79325 (2013).
- ⁶Arndt, S., Unger, P., Berneburg, M., Bosserhoff, A.-K. & Karrer, S. Cold atmospheric plasma (CAP) activates angiogenesis-related molecules in skin keratinocytes, fibroblasts and endothelial cells and improves wound angiogenesis in an autocrine and paracrine mode. *J Dermatol Sci* 89, 181–190 (2018).
- ⁷Boekema, B. et al. Antibacterial and safety tests of a flexible cold atmospheric plasma device for the stimulation of wound healing. *Appl Microbiol Biotechnol* 105, 2057–2070 (2021).
- ⁸Bolgeo, T. et al. The Role of Cold Atmospheric Plasma in Wound Healing Processes in Critically Ill Patients. *J Pers Med* 13, 736 (2023).
- ⁹Xin, J. et al. Effect of cold atmospheric plasma on common oral pathogenic microorganisms: a narrative review. *Ann Med* 57, 2457518.
- ¹⁰Isbary, G. et al. Cold atmospheric plasma devices for medical issues. *Expert Rev Med Devices* 10, 367–377 (2013).
- ¹¹Heinlin, J. et al. Contact-free inactivation of *Trichophyton rubrum* and *Microsporum canis* by cold atmospheric plasma treatment. *Future Microbiol* 8, 1097–1106 (2013).
- ¹²Daeschlein, G. et al. In Vitro Killing of Clinical Fungal Strains by Low-Temperature Atmospheric-Pressure Plasma Jet. *IEEE Transactions on Plasma Science* 39, 815–821 (2011).
- ¹³Braun, J. et al. SARS-CoV-2-reactive T cells in healthy donors and patients with COVID-19. *Nature* <https://doi.org/10.1038/s41586-020-2598-9> (2020) doi:10.1038/s41586-020-2598-9.
- ¹⁴Bekeschus, S., Kramer, A., Suffredini, E., von Woedtke, T. & Colombo, V. Gas Plasma Technology—An Asset to Healthcare During Viral Pandemics Such as the COVID-19 Crisis? *IEEE Transactions on Radiation and Plasma Medical Sciences* 4, 391–399 (2020).
- ¹⁵Arjunan, K. P., Sharma, V. K. & Ptasincka, S. Effects of Atmospheric Pressure Plasmas on Isolated and Cellular DNA—A Review. *International Journal of Molecular Sciences* 16, 2971–3016 (2015).
- ¹⁶Bernhardt, T. et al. Plasma Medicine: Applications of Cold Atmospheric Pressure Plasma in Dermatology. *Oxidative Medicine and Cellular Longevity* 2019, 1–10 (2019).
- ¹⁷Balzer, J. et al. Non-Thermal Dielectric Barrier Discharge (DBD) Effects on Proliferation and Differentiation of Human Fibroblasts Are Primary Mediated by Hydrogen Peroxide. *PLOS ONE* 10, e0144968 (2015).
- ¹⁸Xu, G.-M. et al. Dual effects of atmospheric pressure plasma jet on skin wound healing of mice: Cold plasma promotes wound healing. *Wound Repair and Regeneration* 23, 878–884 (2015).
- ¹⁹Isbary, G. et al. Cold atmospheric argon plasma treatment may accelerate wound healing in chronic wounds: Results of an open retrospective randomized controlled study in vivo. *Clinical Plasma Medicine* 1, 25–30 (2013).
- ²⁰Strohal, R., Dietrich, S., Mittlböck, M. & Hämmerle, G. Chronic wounds treated with cold atmospheric plasmajet versus best practice wound dressings: a multicenter, randomized, non-inferiority trial. *Sci Rep* 12, 3645 (2022). (KRINKO) beim Robert Koch-Institut. *Bundesgesundheitsbl* 61, 448–473 (2018).
- ²¹Hämmerle, G., Ascher, S. & Gebhardt, L. Positive effects of cold atmospheric plasma on pH in wounds: a pilot study. *J Wound Care* 32, 530–536 (2023).
- ²²Strohal, R., Mittlböck, M., Gebhardt, L. & Hämmerle, G. Treatment of chronic wounds with cold plasma: a randomised, single-blind, placebo-controlled clinical study. *J Wound Care* 1–13 (2025) doi:10.12968/jowc.2025.0207.

- ²³ Gnat, S. et al. Cold atmospheric pressure plasma (CAPP) as a new alternative treatment method for onychomycosis caused by *Trichophyton verrucosum*: in vitro studies. *Infection* 49, 1233–1240 (2021).
- ²⁴ Shapourzadeh, A. et al. Inhibitory effects of cold atmospheric plasma on the growth, ergosterol biosynthesis, and keratinase activity in *Trichophyton rubrum*. *Archives of Biochemistry and Biophysics* 608, 27–33 (2016).
- ²⁵ Terabe, Y., Kaneko, N., Nakabayashi, K. & Ando, H. Using Plasma to Treat Chronic Foot Ulcer Infections. *International Journal of Surgical Wound Care* 3, 33–36 (2022).
- ²⁶ Bunz, O. et al. Cold atmospheric plasma as antiviral therapy – effect on human herpes simplex virus type 1. *Journal of General Virology* 101, 208–215 (2020).
- ²⁷ Isbary, G. et al. Randomized placebo-controlled clinical trial showed cold atmospheric argon plasma relieved acute pain and accelerated healing in herpes zoster. *Clinical Plasma Medicine* 2, 50–55 (2014).
- ²⁸ Sutter, J., Bruggeman, P. J., Wigdahl, B., Krebs, F. C. & Miller, V. Manipulation of Oxidative Stress Responses by Non-Thermal Plasma to Treat Herpes Simplex Virus Type 1 Infection and Disease. *Int J Mol Sci* 24, 4673 (2023).
- ²⁹ Sutter, J., Hope, J. L., Wigdahl, B., Miller, V. & Krebs, F. C. Immunological Control of Herpes Simplex Virus Type 1 Infection: A Non-Thermal Plasma-Based Approach. *Viruses* 17, 600 (2025).
- ³⁰ Wang, J. et al. Efficacy and Safety of Cold Atmospheric Plasma-Assisted Therapy for Herpes Zoster: A Randomized, Parallel, Positive-Controlled, Non-inferiority Multicenter Clinical Trial. *Dermatol Ther (Heidelb)* 15, 2391–2408 (2025).
- ³¹ Filipić, A., Gutierrez-Aguirre, I., Primc, G., Mozetič, M. & Dobnik, D. Cold Plasma, a New Hope in the Field of Virus Inactivation. *Trends in Biotechnology* 38, 1278–1291 (2020).
- ³² Bekeschus, S., Schmidt, A., Weltmann, K.-D. & von Woedtke, T. The plasma jet kINPen – A powerful tool for wound healing. *Clinical Plasma Medicine* 4, 19–28 (2016).
- ³³ Neonatal and pediatric wound care A contemporary perspective on innovations and best practices Guido CIPRANDI et. al. - Minerva Medica - Books - Medicine - Pediatrics. <https://www.minervamedica.it/en/books/medical-specialties/pediatrics/scheda.php?cod=L10367>.
- ³⁴ Bae, M., Lademann, J., Meinke, M. C., Meder, B. & Geilen, C. Therapeutic Use of Cold Atmospheric Plasma for the Treatment of Mild Acne Papulopustulosa—A Randomized, Controlled, Double-Blind Pilot Study. *Dermatologic Therapy* 2025, 4228323 (2025).
- ³⁵ Assadian, O. et al. Effects and safety of atmospheric low-temperature plasma on bacterial reduction in chronic wounds and wound size reduction: A systematic review and meta-analysis. *Int Wound J* 16, 103–111 (2019).
- ³⁶ Mueller, R. S., Bergvall, K., Besignor, E. & Bond, R. A review of topical therapy for skin infections with bacteria and yeast. *Vet Dermatol* 23, 330–341, e62 (2012).
- ³⁷ Emmert, S. Plasmamedizin - eine Innovation weit über die Dermatologie hinaus. *JDDG: Journal der Deutschen Dermatologischen Gesellschaft* 13, 95–96 (2015).
- ³⁸ Boeckmann, L. et al. Aktuelle Indikationen der Plasmatherapie in der Dermatologie. *Hautarzt* 71, 109–113 (2020).
- ³⁹ Karrer, S. et al. Atmospheric Plasma in the Treatment of Acne Vulgaris. *Applied Sciences* 11, (2021).
- ⁴⁰ Dejonckheere, C. S. et al. Non-Invasive Physical Plasma for Preventing Radiation Dermatitis in Breast Cancer: A First-In-Human Feasibility Study. *Pharmaceutics* 14, 1767 (2022).
- ⁴¹ Dejonckheere, C. S. et al. Non-invasive physical plasma for preventing radiation dermatitis in breast cancer: Results from an inpatient-randomised double-blind placebo-controlled trial. *Clinical and Translational Radiation Oncology* 44, 100699 (2024).
- ⁴² Dejonckheere, C. S., Schmeel, L. C. & Stope, M. B. Concise review of non-invasive physical plasma as a promising treatment option for radiation injuries of the skin. *Wound Repair Regen* 31, 415–417 (2023).
- ⁴³ Bai, F., Ran, Y., Zhai, S. & Xia, Y. Cold Atmospheric Plasma: A Promising and Safe Therapeutic Strategy for Atopic Dermatitis. *International Archives of Allergy and Immunology* 1–14 (2023) doi:10.1159/000531967.

⁴⁴ Arisi, M. et al. Cold Atmospheric Plasma (CAP) for the Treatment of Actinic Keratosis and Skin Field Cancerization: Clinical and High-Frequency Ultrasound Evaluation. *Dermatol Ther (Heidelb)* <https://doi.org/10.1007/s13555-021-00514-y> (2021) doi:10.1007/s13555-021-00514-y.

⁴⁵ Friedman, P. C. Cold atmospheric pressure (physical) plasma in dermatology: where are we today? *International Journal of Dermatology* 59, 1171–1184 (2020).

⁴⁶ Friedman, P. C., Miller, V., Fridman, G. & Fridman, A. Use of cold atmospheric pressure plasma to treat warts: a potential therapeutic option. *Clinical and Experimental Dermatology* 44, 459–461 (2019).

⁴⁷ Friedman, P. C., Miller, V., Fridman, G., Lin, A. & Fridman, A. Successful treatment of actinic keratoses using nonthermal atmospheric pressure plasma: A case series. *Journal of the American Academy of Dermatology* 76, 349–350 (2017).

⁴⁸ Friedman, P. C., Miller, V., Fridman, G. & Fridman, A. Use of cold atmospheric pressure plasma to treat warts: a potential therapeutic option. *Clin Exp Dermatol* 44, 459–461 (2019).

⁴⁹ Friedman, P. C., Fridman, G. & Fridman, A. Cold atmospheric plasma device induces hair growth: A nonrandomized, controlled study. *J Am Acad Dermatol* 93, 294–295 (2025).

⁵⁰ Koch, F. et al. Efficacy of cold atmospheric plasma vs. diclofenac 3% gel in patients with actinic keratoses: a prospective, randomized and blindfolded study (ACTICAP). *Journal of the European Academy of Dermatology and Venereology* <https://doi.org/10.1111/jdv.16735> (2020) doi:10.1111/jdv.16735.

⁵¹ Terabe, Y., Kaneko, N. & Ando, H. Treating hard-to-heal skin and nail onychomycosis of diabetic foot with plasma therapy. *Dermatol Ther* 34, e15127 (2021).

⁵² Misra, N. N., Yadav, B., Roopesh, M. S. & Jo, C. Cold Plasma for Effective Fungal and Mycotoxin Control in Foods: Mechanisms, Inactivation Effects, and Applications: Cold plasma for effective fungal.... *Comprehensive Reviews in Food Science and Food Safety* 18, 106–120 (2019).

⁵³ Preissner, S. et al. Adjuvant antifungal therapy using tissue tolerable plasma on oral mucosa and removable dentures in oral candidiasis patients: a randomised double-blinded split-mouth pilot study. *Mycoses* 59, 467–475 (2016).

⁵⁴ Roberts, D. M. et al. Cold atmospheric plasma improves antifungal responsiveness of *Aspergillus flavus* and *Fusarium keratoplasticum* conidia and mycelia. *PLoS One* 20, e0326940 (2025).

⁵⁵ Sun, P. P. et al. Atmospheric pressure cold plasma as an antifungal therapy. *Applied Physics Letters* 98, 021501 (2011).

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terraplasma medical GmbH