



DERMA  
plasma care<sup>®</sup>

Cold plasma therapy

CAP for treatment of bacterial/fungal skin infections



SAVE AND  
RELIABLE



SUITABLE FOR PATIENTS  
WITH PACEMAKERS



PORTABLE AND  
EASY-TO-USE



**NAIL MYCOSIS**  
initial situation &  
end of treatment




**BURNS**  
initial situation & after  
5 CAP treatments



**FINGERTIP NECROSIS**  
initial situation &  
end of treatment

**ACTINIC KERATOSIS**  
initial situation &  
end of treatment

**LEG ULCER**  
initial situation &  
end of treatment



## plasma derma care – new chances in dermatology

**The plasma derma care is a mobile medical device for the treatment of acute and chronic bacterial or fungal skin diseases or certain skin dysplasias/neoplasias such as actinic keratosis.**

This portable device generates cold atmospheric plasma (CAP) by supplying high energy to the ambient air. CAP effectively inactivates bacteria (independent of type and level of antibiotic resistance) and human pathogenic fungi.

Furthermore, CAP has other positive effects that promote healing of dermatological diseases – e.g. an antipruritic (anti-itch)<sup>1</sup> effect. By stimulating the immune system locally, CAP is potentially effective for treatment of precancerous skin, too.<sup>2,3</sup>

Most dermatological diseases are caused by pathogens such as bacteria and fungi or are associated with those. Patients often suffer from a high level of psychological strain. Apart from this, various pathogens can also lead to a more severe course of disease.

Thus, early and effective treatment of bacterial/fungal skin diseases is extremely important for a successful therapy.

The use of **plasma derma care** provides a completely new dermatological therapy approach for the effective treatment of skin diseases without side effects.

Cover picture shows female patient ▲  
with severe atopic dermatitis on both hands and wrists.

# Our vision of cold plasma in medicine

**Terraplasma medical GmbH looks back on many years of experience in the development of plasma technologies and their application in various fields.**

As a spin-off of Max Planck Institute for Extraterrestrial Physics headed by Prof. Dr. Dr. Gregor Morfill and Dr. Julia Zimmermann, **terrapiasma GmbH** was founded in 2011.

**Terrapiasma GmbH** still focuses on the development of various plasma technologies for applications in medical technology, hygiene and air pollution control.

In 2016, **terrapiasma medical GmbH** was founded as an independent subsidiary which since has developed and researched medical devices based on plasma technology. The **plasma derma care**® uses cold atmospheric plasma designed by terrapiasma medical GmbH and already applied in various medical devices to inactivate microorganisms (bacteria incl. MRE and human pathogenic fungi).

For years, plasma technology has already been established in wound treatment through the use of **plasma care**®.

Apart from this, another medical device – the **plasma intensive care**® – was pre-launched at the beginning of 2021 for treatment of ventilated COVID-19 patients.

**The main goal of terrapiasma medical GmbH is to convince the medical world of the innovative and promising plasma technology and to support patients in activating the body's own healing processes through a modern therapeutic approach - and this against the background of increasing antibiotic resistances and allergies.**



# Effects of cold plasma on bacteria and viruses

**Cold atmospheric plasma (CAP) is known to have a strong antimicrobial and antiviral effect while preserving human tissues and cells.**

The active components (reactive species) generated by CAP can interact with cells in various ways. The effects are physical (e. g. recombination or de-excitation of excited molecules/atoms at the surface) as well as chemical (e. g. hydrogen denaturation by reactions of hydroxyl radicals).

In bacteria (prokaryotic cells) – or more precisely their cellular macro-molecules (including DNA) – a direct destruction takes place.

**This effect also applies for bacteria with antibiotic resistance.**<sup>4,5,6</sup> Furthermore, an antiviral effect of the reactive oxygen species produced by cold atmospheric plasma has been observed.<sup>7</sup> This is presumably based on a modification of viral proteins leading to the inactivation of the viruses. In eukaryotic cells, such as human cells, DNA is protected by nucleus and its membrane as well as cyto-logical repair mechanisms. In addition, those cells are protected within cell compounds. This means that there is no risk of cell damage.


Quite the opposite, **in human cells a local immune response is triggered and cell growth and division are stimulated**, thus promoting the regeneration of skin or mucous membrane.<sup>8,9,10</sup>

This effect is locally limited to the uppermost cell layer. Among other reasons, this is due to a very short half-life of the reactive species, which will have already completely reacted by interacting with the upper cells or surrounding cell milieu.<sup>11,12</sup>

**Cold atmospheric plasma (CAP) does not cause any un-specific cell damage due to its low temperature and therefore is generally suitable for medical application.**

**1. INACTIVATION OF BACTERIA + VIRUSES**  
**CAP interacts with microorganisms** (e.g. bacteria + viruses), inactivates the cell and prevents cell proliferation by destroying the DNA.

**2. ACTIVATION OF (WOUND) HEALING**  
**Cell separation** in healthy cells **is stimulated by CAP**. Subsequent intracellular processes: release of cytokines, activation of angiogenesis and metabolism.



▲ Figure 3: effects of cold plasma



# Effects of cold plasma on fungi

Fungal infections are also a common cause for skin diseases. Some types of fungi, such as *Candida albicans*, are already naturally present on the skin, but can cause an infection, if the immune system is weakened.


Apart from inactivating bacteria (including multi-drug-resistant pathogens)<sup>13,14,15,16</sup>, CAP can also kill fungi.<sup>17,18,19</sup> Like human cells, fungi belong to the eukaryotes, in contrast to bacteria. They thus possess a nucleus in which the DNA is located.

Nevertheless, fungi can be killed by reactive oxygen species (ROS) – a mechanism that is also used by the immune system to fight *Candida albicans* infections (fig.4).<sup>20</sup>

Treatment with CAP causes deformation of fungal spores, which leads to their rupture, flattening and shrinkage. Also, the DNA within the spores can be destroyed. Reactive oxygen species (ROS) trigger various reactions in fungi. Depending on the dose, these can lead to the oxidation of intracellular membranes and proteins, to structural changes within the cell, and finally to apoptosis (programmed cell death).<sup>21</sup>



▲ Figure 4: *Candida albicans* under a microscope



▲ Figure 5: effects of CAP on fungi

## Recent study results confirm high effectiveness

Preclinical evaluations prove the antimicrobial and anti-fungal effectiveness of the **plasma derma care**<sup>®</sup> (see table).

Most of the tested microorganisms are resident on our skin and thus are predestined to cause certain infections. The development of perioral dermatitis, for example, is mostly triggered by *fusobacteriales*.<sup>22,23</sup> Acne, on the other hand, is mainly caused by *propionibacterium acnes*.<sup>24</sup> *Staphylococcus aureus* is a highly resistant bacterial strain that causes nosocomial infections<sup>25</sup> and is often identifiable in bacterial superinfections in atopic dermatitis as well as many other skin diseases.<sup>26</sup>

The following table (table 1) clearly reflects the high effectiveness of CAP produced by **plasma derma care** device against human pathogenic bacteria and fungi. On agar, 99.999 % of tested bacteria and of the yeast *C. albicans* were killed within 60 seconds.

Moreover, the **plasma derma care**<sup>®</sup> was effective upon application to biofilms of *E. faecalis* (99.9 % reduction within 1 minute). Under more life-like conditions in ex vivo models of pig skin and human skin, 69 to 83 % of bacteria were killed within the same period.

Apart from the described studies of effectiveness, **plasma derma care**<sup>®</sup> has also been thoroughly tested with regard to its safe use and the occurrence of side effects. However, plasma treatments of up to 3 minutes had no impact on vitality, viability or migration behavior of primary fibroblasts and keratinocytes. Nor did "normal" or "sensitive" skin from healthy donor biopsies display any histological or pro-apoptotic changes. Mutagenicity studies (HGPRT test using V79 cells) provided no evidence of any genotoxic potential of CAP, too.

The results of these preclinical studies by **terraplasma group** also correspond to the study results of recent publications.

	pathogenes	achieved log reduction after 1 min plasma treatment	associated skin diseases
bacteria	Fusobacterium Nucleatum	4,4	perioral dermatitis
	s. aureus / MRSA	5,6 / 6,0	acne, neurodermatitis
	p. acnes	.....	acne
	p. aeruginosa	5,9	acne
	s. agalactiae	6,5	Various superficial skin infections, wound infections, bacterial superinfections
	s. constellatus	5,8	
	e. coli	5,8	
	s. constellatus	5,7	
fungi	c.albicans	4,6	candidiasis
	t. rubrum	.....	dermatophytoses
	m. furfur	.....	tinea versicolor

▲ Table 1: Investigated human pathogens with achieved log reduction

**Acne and actinic keratosis** have already successfully been treated with cold atmospheric plasma.<sup>27,28</sup> CAP therapy was found to be an extraordinary effective treatment for patients with actinic keratosis – a non-bacterial/fungal associated skin disease. This effect is presumably based on the activation of local immune processes and the destruction of degenerated cells by the plasma – as already mentioned on left page.

Reports about tissue sensitivity or side effects caused by plasma treatment have not yet occurred in studies.<sup>29</sup> Moreover, the **effectiveness of plasma treatment in bacteria appears to be independent from their level of resistance to antibiotics or the species type**.<sup>30,31,32</sup>

## A wide range of indications: from acne to atopic dermatitis

The multiple mechanism of plasma action provides treatment for a wide range of dermatological diseases (see table 2). Due to its flexible biocompatible foam layer, even hard-to-reach treatment areas could be easily accessed by plasma derma care® spacer.

**Mycoses** and **bacteriosis** are some of the most common dermatological diseases and – in case of delayed or inadequate therapy – they may develop severe clinical symptoms.

It is also worth to emphasize that **acne vulgaris**, with a prevalence of about 70-95%, is among adolescents<sup>33</sup> the most frequent dermatological disease worldwide.

Currently, **bacterial infections** are mainly treated with antibiotics. However, prolonged antibiotic treatment is not recommended due to the risk of developing resistances. Mycoses are usually treated with topical antifungal agents. Therapy may cause side effects such as skin irritation, redness or hypersensitivity reactions; recurrences are very common.

Due to its antibacterial and antifungal effect, treatment with **plasma derma care®** is a suitable alternative for the local therapy of mycoses and bacteriosis.

The treatment of **atopic dermatitis** as a chronic constitutional disease is difficult. The usual therapies with glucocorticoids or antibiotics should not be continued for a longer time. In contrast, treatment with cold plasma is also possible in the long term, since no development of resistances or further side effects are to be expected.

**Actinic keratosis** can be treated well in early stages with current treatment methods such as photodynamic therapy (PDT). However these methods are often associated with recurrences and side effects such as severe pain symptoms<sup>34</sup>. These treatment methods are often used in combination to have synergistic effects in combating the disease. At this point, **plasma derma care®** can be used as an add-on therapy.

	INDICATIONS	ESTABLISHED THERAPIES	plasma derma care
<b>bacterial infections (bacteriosis)</b>	akne vulgaris	antibiotics, antiseptics, glucocorticoids, retinoids, azelaic acid & benzoyl peroxide (for acne vulgaris)	antibacterial, promotes wound healing, prevents superinfection and progression of inflammatory reactions ⊕ no resistances ⊕ long-term use possible
	perioral dermatitis		
	soft-tissue infections (folliculitis, boils, abscesses)	⊖ retinoids with strong teratogenic effect ⊖ antibiotics involve the risk of resistances ⊖ glucocorticoids can only be used for short periods	
	impetigo contagiosa		
<b>fungal infections (mycosis)</b>	pityriasis versicolor (tinea versicolor)	antifungal agents (azoles, polyenes, allylamines, etc.) – when used locally often associated with skin redness and irritations	antimycotic alternative/ combination therapies ⊕ prevention of recurrences and superinfection
	dermatophytosis (tinea)		
<b>other skin diseases</b>	atopic dermatitis (neurodermatitis)	calcineurin antagonists, glucocorticoids, antibiotics, hyposensitization	antipruritic, antibacterial, wound healing (in case of excoriations (skin-picking))
	actinic keratosis	cryotherapy, 5-fluorouracil, laser therapy, photodynamic therapy (PDT) ⊖ severe pain symptoms especially with PDT	activation of local immune processes, destruction of degenerated cells ⊕ without side effects, treatment painless
	viral skin diseases with bacterial superinfection (e.g. herpes labiales)	antivirals (aciclovir, ganciclovir, etc.), antibiotics (for bacterial superinfections)	in combination therapy: for the prevention/ treatment of bacterial superinfections (antipruritic and analgetic effect)

▲ Table 2: Indications for plasma derma care

## How to use – simple, effective and painless\*

The **plasma derma care** is a portable, battery-operated and therefore mobile medical device. Accordingly, it can easily be used in the clinical or practice sector, but also by outpatient care services. Treatment with the **plasma derma care** only requires a few minutes and its use is very intuitive.

Plasma therapy does not exclude treatment with other therapeutic agents. Topical/systemic glucocorticoids, retinoids, phototherapy, etc. remain an effective add-on and may also be used depending on the individual patient profile. However, topically applied preparations should be removed prior to plasma treatment so that CAP can freely react.

In general, the use of the device is similar for all indications listed in table 2.


The affected skin area is cleaned prior to treatment and any existing crusts or pus are removed before the flexible biocompatible foam spacer is placed on the skin. The gaseous plasma can perfectly diffuse to the affected areas and unfold its effects on the skin surface (see figure 4).

An antiseptic can then be applied to provide a long-term antibacterial effect. The recommended treatment time is 1 -3 minutes per affected area.

The affected skin area is only exposed to minimal contact (application of the spacer), which means a reduced pain potential compared to other topical treatment options. Areas close to eyes and mouth/nose openings should be treated very cautiously.

Plasma should not enter the respiratory tract or get into direct contact with the eyes; if necessary, aids such as nasal plugs and protective glasses can be used.

*\*Patients described treatment with the **plasma derma care** as free of pain.*



▲ Figure 4: treatment process with the **plasma derma care**



## The advantages at a glance

Treatment with **plasma derma care** represents a completely new therapeutic approach for diverse chronic and acute skin diseases and offers a lot of advantages in comparison to conventional therapeutics.

In particular, there are no contraindications, the application is free of side effects, and at the same time treatment is highly effective.

Due to the 2-fold effect, i. e. activation of a local immune response and antimicrobial as well as antipruriginous effects, CAP can contribute in different ways to symptomatic relief or recovery.

Treatment with **plasma derma care** is quick and easy to perform, painless, cost effective and economical. Thus it may effectively be added to or even replace conventional therapies.

- high antiseptic efficacy
- wide range of indications for dermatology – from acne to neurodermatitis
- broad spectrum of effects – inactivation of all bacteria on skin
- no adverse events, skin intolerances and long-term risks <sup>29</sup>
- no resistances <sup>6</sup> or allergies
- treatment with minimum wound contact (spacer) and thus considerable reduction of pain potential in comparison to topical antiseptics.
- CAP inactivates bacteria - independent from their resistance to antibiotics or their species <sup>30,31,32</sup>
- also hard-to-reach or facial skin can be easily treated with flexible foam spacer

## Case report

### WOUND HEALING DISORDER – BASAL CELL CARCINOMA

**82-year-old male patient suffers from basal cell carcinoma on parietal bone:** removal of basal cell carcinoma in a doctor's practice. Delay in wound healing and occurrence of first signs of infection. Formation of necrosis at front wound edge. Wound with pain sensitivity to pressure.

#### Treatment:

**Picture 1:** Start of CAP therapy on April 28th, 2022 in order to avoid further complications. After debridement with a physiological saline solution wound was treated 2 min with CAP 2 times per week.

**Picture 2:** Visible improvement of wound healing during 2nd treatment week. Wound edges no longer necrotic and no more visible signs of infection. Patient is free of pain.

**Picture 3:** Complete epithelization of wound on June 14<sup>th</sup>, 2022 after 11 CAP treatments within 6 weeks.



Initial situation: first CAP treatment



Course of treatment



After 11 CAP treatments

# Case report

## SYMPTOMATIC TREATMENT OF PSORIASIS

### 52-year-old male patient suffering from psoriasis for 6 months.

Psoriasis was initially treated with cortisone ointment with short-term improvement of symptoms. The skin soon became itchy, inflammatory and scaly again.

Skin care products used: urea ointment 10% and salicylic acid patch.

#### Treatment:

Picture 1: Heavy symptoms with bleeding contact dermatitis at right ankle. Skin areas were treated with 1 min of CAP-therapy each.

Picture 2: Patient reported a clear, subjective improvement of symptoms. The itching was reduced and the strained feeling of skin decreased. The open, scratched skin areas were found to be almost completely closed. Skin redness was visibly reduced.

Picture 3: After only 1 week and 4 KAP treatments, the patient's subjective sensation turned out to be very positive: no pain, hardly any itching, no more feeling of skin tension on the ankle. Objectively, there was a significant reduction of inflammation, skin condition improved and contact dermatitis was almost gone.



Initial situation: 1<sup>st</sup> CAP treatment



Course of treatment



After one week/ 4 CAP treatments

# Case report

## AKNE INVERSA

### 44-year-old female patient suffering from acne inversa.

Severe wound defect in patient with acne inversa around right axilla. Large wound areas up to 4.5 cm in diameter with multiple fistula tracts and undermined wound edges. Heavily exuding wounds on scar tissue due to multiple recurrences.

#### Treatment:

15 CAP treatments within 8 weeks. The complete wound area was divided into sections. Each section was treated with CAP for 1 minute each (6 x 1 minute). Wounds dried out significantly after 10 CAP treatments and wound edges smoothed.

After 15 CAP treatments, there was hardly any wound exudate and no more wound edge undermining. The wounds closed gradually.

The therapeutic intervals were reduced to 1 treatment per week. A final result is still pending.



Initial situation: 1<sup>st</sup> CAP treatment



After 10 CAP treatments



After 15 CAP treatments

# Case report

## EXTRACT FROM "DERMATOLOGIC THERAPY", WILEY, AUGUST 27<sup>TH</sup>, 2021

### 80-year old female patient suffering from diabetes mellitus (HbA1c 8,5%)

Past medical history: amputation of left foot below the knee and amputation of the right fourth toe, 3 years ago.

The patient developed an onychomycosis on the skin and toenails of right foot.

After 4 weeks of non-successful antifungal and antibiotic treatment, the patient was treated with CAP for 10 min on a daily basis. After only one week a clear improvement could be observed.

2 weeks later topical therapy was changed from an antifungal ointment to a moisturizing ointment with supplemental plasma therapy (fig. B).

Patient finished plasma therapy after 3 weeks in total (fig.C).



Read the whole article:



(A) Before treatment



(B) After two weeks of CAP treatment




(C) After three weeks of CAP treatment

Those who suffer from nail mycosis are often ashamed and hide their nails in shoes and socks. Medically, a therapy of nail mycosis is advisable to prevent a bacterial superinfection (e.g. erysipelas) and a spreading skin mycosis.

### In case of diabetic foot syndrome early nail mycosis therapy is indispensable!

In more severe cases, patients must also receive a systemic antimycotic therapy. The side effects of this therapy are often described by patients as very uncomfortable. Treatment with CAP can not only eliminate the fungal disease, but also prevent its spread and reinfection.

The advantages in comparison to conventional therapies are quite obvious.




# Excursus: NAIL MYCOSIS

## First reports about treatment of nail mycosis


### Frequency of treatment:

Every toenail is different: basically, a nail should be retreated until a completely healthy nail has been re-grown. In previous case reports, patients have been treated at least 1 to 2 times per week.

Once the healthy nail progresses to re-grow, the frequency of therapy sessions can be reduced while monitoring.



Initial situation



4 weeks/ 6 CAP treatments



Initial situation



End of CAP therapy

# Literaturverzeichnis

- 1 Heinlin, J. et al. A randomized two-sided placebo-controlled study on the efficacy and safety of atmospheric non thermal argon plasma for pruritus. *J Eur Acad Dermatol Venereol* 27, 324–331 (2013).
- 2 Friedman, P. C., Miller, V., Fridman, G. & Fridman, A. Use of cold atmospheric pressure plasma to treat warts: a potential therapeutic option. *Clin Exp Dermatol* 44, 459–461 (2019).
- 3 Wirtz, M., Stoffels, I., Dissemond, J., Schadendorf, D. & Roesch, A. Actinic keratoses treated with cold atmospheric plasma. *Journal of the European Academy of Dermatology and Venereology* 32, e37–e39 (2018).
- 4 Nicol, M. J. et al. Antibacterial effects of low-temperature plasma generated by atmospheric-pressure plasma jet are mediated by reactive oxygen species. *Scientific Reports* 10, 3066 (2020).
- 5 Hoon Park, J. et al. A comparative study for the inactivation of multidrug resistance bacteria using dielectric barrier discharge and nano-second pulsed plasma. *Scientific Reports* 5, 13849 (2015).
- 6 Zimmermann, J. L. et al. Test for bacterial resistance build-up against plasma treatment. *New Journal of Physics* 14, 073037 (2012).
- 7 Filipić, A., Gutierrez-Aguirre, I., Primc, G., Mozetič, M. & Dobnik, D. Cold Plasma, a New Hope in the Field of Virus Inactivation. *Trends in Biotechnology* 38, 1278–1291 (2020).
- 8 Arndt, S. et al. Effects of cold atmospheric plasma (CAP) on  $\beta$ -defensins, inflammatory cytokines, and apoptosis-related molecules in keratinocytes in vitro and in vivo. *PLoS ONE* 10, e0120041 (2015).
- 9 Arndt, S. et al. Cold atmospheric plasma (CAP) changes gene expression of key molecules of the wound healing machinery and improves wound healing in vitro and in vivo. *PLoS ONE* 8, e79325 (2013).
- 10 Schmidt, A., Bekeschus, S., Wende, K., Vollmar, B. & von Woedtke, T. A cold plasma jet accelerates wound healing in a murine model of full-thickness skin wounds. *Experimental Dermatology* 26, 156–162 (2017).
- 11 Leduc, M., Guay, D., Leask, R. L. & Coulombe, S. Cell permeabilization using a non-thermal plasma. *New Journal of Physics* 11, 115021 (2009).
- 12 Von Woedtke, T., Schmidt, A., Bekeschus, S., Wende, K. & Weltmann, K.-D. Plasma Medicine: A Field of Applied Redox Biology. *In Vivo* 33, 1011–1026 (2019).
- 13 Isbary, G.; Morfill, G.; Schmidt, H. U.; Georgi, M.; Ramrath, K.; Heinlin, J.; Karrer, S.; Landthaler, M.; Shimizu, T.; Steffes, B.; Bunk, W.; Monetti, R.;
- 14 Zimmermann, J. L.; Shimizu, T.; Schmidt, H.-U.; Li, Y.-F.; Morfill, G. E.; Isbary, G. Test for Bacterial Resistance Build-up against Plasma Treatment. *New J. Phys.* 2012, 14 (7), 073037. <https://doi.org/10.1088/1367-2630/14/7/073037>.
- 15 Maisch, T.; Shimizu, T.; Li, Y.-F.; Heinlin, J.; Karrer, S.; Morfill, G.; Zimmermann, J. L. Decolonisation of MRSA, S. Aureus and E. Coli by Cold-Atmospheric Plasma Using a Porcine Skin Model in Vitro. *PLoS One* 2012, 7 (4), e34610. <https://doi.org/10.1371/journal.pone.0034610>.
- 16 Becker, S.; Zimmermann, J. L.; Baumeister, P.; Brunner, T. F.; Shimizu, T.; Li, Y.-F.; Morfill, G. E.; Harréus, U.; Welz, C. Effects of Cold Atmospheric Plasma (CAP) on Bacteria and Mucosa of the Upper Aerodigestive Tract. *Auris. Nasus. Larynx* 2019, 46 (2), 294–301. <https://doi.org/10.1016/j.anl.2018.07.008>.
- 17 Heinlin, J.; Maisch, T.; Zimmermann, J. L.; Shimizu, T.; Holzmann, T.; Simon, M.; Heider, J.; Landthaler, M.; Morfill, G.; Karrer, S. Contact-Free Inactivation of *Trichophyton Rubrum* and *Microsporum Canis* by Cold Atmospheric Plasma Treatment. *Future Microbiol.* 2013, 8 (9), 1097–1106. <https://doi.org/10.2217/fmb.13.86>.
- 18 Shapourzadeh, A.; Rahimi-Verki, N.; Atyabi, S.-M.; Shams-Ghahfarokhi, M.; Jahanshiri, Z.; Irani, S.; Razzaghi-Abyaneh, M. Inhibitory Effects of Cold Atmospheric Plasma on the Growth, Ergosterol Biosynthesis, and Keratinase Activity in *Trichophyton Rubrum*. *Arch. Biochem. Biophys.* 2016, 608, 27–33. <https://doi.org/10.1016/j.abb.2016.07.012>.
- 19 Daeschlein, G.; Scholz, S.; von Woedtke, T.; Niggemeier, M.; Kindel, E.; Brandenburg, R.; Weltmann, K.-D.; Junger, M. In Vitro Killing of Clinical Fungal Strains by Low-Temperature Atmospheric-Pressure Plasma Jet. *IEEE Trans. Plasma Sci.* 2011, 39 (2), 815–821. <https://doi.org/10.1109/TPS.2010.2063441>.
- 20 Kaloriti, D.; Jacobsen, M.; Yin, Z.; Patterson, M.; Tillmann, A.; Smith, D. A.; Cook, E.; You, T.; Grimm, M. J.; Bohovych, I.; Grebogi, C.; Segal, B. H.; Gow, N. A. R.; Haynes, K.; Quinn, J.; Brown, A. J. P. Mechanisms Underlying the Exquisite Sensitivity of *Candida Albicans* to Combinatorial Cationic and Oxidative Stress That Enhances the Potent Fungicidal Activity of Phagocytes. *mBio* 2014. <https://doi.org/10.1128/mBio.01334-14>.
- 21 Misra, N. N.; Yadav, B.; Roopesh, M. S.; Jo, C. Cold Plasma for Effective Fungal and Mycotoxin Control in Foods: Mechanisms, Inactivation Effects, and Applications: Cold Plasma for Effective Fungal.... *Compr. Rev. Food Sci. Food Saf.* 2019, 18 (1), 106–120. <https://doi.org/10.1111/1541-4337.12398>.
- 22 Takiwaki, H., Tsuda, H., Arase, S. & Takeichi, H. Differences between intrafollicular microorganism profiles in perioral and seborrhoeic dermatitis: Intra-follicular microorganisms in dermatitis. *Clinical and Experimental Dermatology* 28, 531–534 (2003).
- 23 Berardi, P., Benvenuti, S., Genga, A. & Cecchini, F. Demonstration of fusobacteria in eruptions of perioral dermatitis using the tape stripping toluidine blue (TSTB) method. *Journal of the European Academy of Dermatology and Venereology* 3, 495–499 (2006).
- 24 *Dermatologie*. b-002-11376 (Georg Thieme Verlag, 2010). doi:10.1055/b-002-11376.
- 25 Chudobova, D. et al. Influence of microbiome species in hard-to-heal wounds on disease severity and treatment duration. *The Brazilian Journal of Infectious Diseases* 19, 604–613 (2015).
- 26 Goebeler, M. et al. *Basiswissen Dermatologie*. (Springer, 2017).
- 27 Mariachiara, Dr. A. et al. Cold atmospheric plasma (CAP) as a promising therapeutic option for mild to moderate acne vulgaris: Clinical and non-invasive evaluation of two cases. *Clinical Plasma Medicine* 100110 (2020) doi: 10.1016/j.cpme.2020.100110.
- 28 Prof. Dr. med. Rösch, A., Salva, K. & McGovern, M. DRKS00012714 - ACTICAP; Eine prospektive, randomisierte, monozentrische, untersucherblindete Studie zur Evaluation der Veränderung des kutanen Mikrobioms in Korrelation zum Ansprechen von kaltem atmosphärischem Plasma im Vergleich zu Diclofenac 3 % in 2,5 % Hyaluronsäure (Solaraze 3% Gel(R)) in Patienten mit aktinischen Keratosen. [https://www.drks.de/drks\\_web/navigate.do?navigationId=trial.HTML&TRIAL\\_ID=DRKS00012714](https://www.drks.de/drks_web/navigate.do?navigationId=trial.HTML&TRIAL_ID=DRKS00012714).
- 29 Kramer, A. et al. Consensus on Wound Antisepsis: Update 2018. *Skin Pharmacology and Physiology* 31, 28–58 (2018).
- 30 Zimmermann, J. L. et al. Test for bacterial resistance build-up against plasma treatment. *New Journal of Physics* 14, 073037 (2012).
- 31 Heinlin, J. et al. Plasma applications in medicine with a special focus on dermatology: Plasma medicine. *Journal of the European Academy of Dermatology and Venereology* 25, 1–11 (2011).
- 32 Maisch, T. et al. Contact-free cold atmospheric plasma treatment of *Deinococcus radiodurans*. *J. Ind. Microbiol. Biotechnol.* 39, 1367–1375 (2012).
- 33 WMF & Deutsche Dermatologische Gesellschaft. S2 Leitlinie - Behandlung der Akne. (2011).
- 34 Fink, C., Enk, A. & Gholam, P. Aspekte der Schmerztherapie bei der photodynamischen Therapie. *JDDG: Journal der Deutschen Dermatologischen Gesellschaft* 13, 15–22 (2015).





**terraplasma medical GmbH**

Parking 32

85748 Garching | Germany

**Phone** + 49 89 588 055 30

**eFax** + 49 89 588 055 99

**Let's stay in touch:**

**Linked in** [terraplasma-medical](https://www.linkedin.com/company/terraplasma-medical)


**Mail** [info@terraplasma-medical.com](mailto:info@terraplasma-medical.com)

**Web** [www.terraplasma-medical.com](http://www.terraplasma-medical.com)

**Jens Kirsch** | Chief Executive Officer (CEO)

**PD Dr. rer. nat. habil. med. Julia Zimmermann** | Chief Medical Officer (CMO)

**Prof. Dr. Dr. h.c. Gregor Morfill** | Scientific advisory board



Stempel



USER-FRIENDLY



NO DEVELOPMENT  
OF RESISTANCES



REDUCTION OF  
BACTERIAL LOAD